

## VWR 1415 Vacuum Oven



VWR 1415 Signature microprocessor-controlled vacuum oven, small capacity, 0.6 cu.ft. Demo model barely used.

Rating: Not Rated Yet

**Price**

Sales price \$1,350.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

### Description

**!!!! - RECENTLY SOLD - !!!!**

VWR 1415 Signature microprocessor-controlled vacuum oven, small capacity, 0.6 cu.ft. LED digital control and readout, built-in 0-30" Hg vacuum gauge. Floating tempered glass viewing window assures a tight seal. Stainless steel chamber. Demo model barely used.