

## VWR 2.3 cu. ft. Vacuum Oven



65.1L (2.3 cu. ft.); 120V; 1600w; 13.3A; LED display

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Manufacturer [VWR](#)

### Description

Model 6292

Capacity: 2.3 cu.ft. / 65.1 L

Chamber Material: Stainless Steel

Temperature (Metric): Ambient +10° to 220°C

Type: OTP

Display: LED with 1°C increments

Material: Powder-coated Heavy-gauge Steel

No. of Shelves: 2

Shelf Material: Aluminum

Temperature Resolution:  $\pm 1.0^{\circ}\text{C}$

Temperature Uniformity:  $\pm 5.0^{\circ}\text{C}$  at 100°C 25 ft.Hg

Vacuum Range: 0 to 30 in.Hg

Dimensions (D x W x H) Interior: 20 x 14 x 14 in. (51 x 36 x 36 cm)

Dimensions (L x W x H) Exterior: 25 x 25 x 22 in. (64 x 64 x 56 cm)

Shipping Weight (English): 295 lb./ 134 kg

Electrical Requirements : 120 V 50/60 Hz

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)

