

## VWR HCAFS1420 General Purpose Auto Defrost Laboratory Freezer, -20





14 cu. ft.

Rating: Not Rated Yet

**Price**

Sales price \$650.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

Description

Auto-Defrost Freezers are ideal for industrial, clinical, or research laboratory requirements. Designed to provide dependable, energy-efficient performance, and manufactured with CFC-free high-density urethane foam cabinet and door insulation to maintain temperature stability. Auto-Defrost design provides the convenience of a frost-free interior. However, chamber air temperatures do experience a moderate rise during defrost cycles.

Specifications:

Capacity: 14 cu. ft.

Exterior dimensions: 28"W x 26 1/4"D x 67 3/4"H

Door: One solid swing door, right hinged

Shelves: Four adjustable shelves and 1 wire basket

Electrical: 5 Amps, NEMA 5-15 Plug, 115VAC, 60Hz

Refrigeration System: R600a

Adjustable operating temperature range: -15°C to -20°C

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)