

## VWR Incubating Microplate Shaker



Temperature Range: Ambient +5...65 °C (149 °F)

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Manufacturer [VWR](#)

### Description

Ideal for cell cultures, hybridizations, and suspensions.

### Specifications

Temperature Range: Ambient +5...65 °C (149 °F)

Temperature Uniformity:  $\pm 0.5^{\circ}\text{C}$  at  $37^{\circ}\text{C}$

Speed accuracy:  $\pm 2\%$

Max. Load: 2.3 kg (5 lbs.)

Overall Dimensions: 28W x 19.7H x 43.2D cm (11 x 73/4 x 17")

Timer: 1 second to 160 hours

Shipping Weight: 13.6 kg (30 lbs.)

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
**[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)**