

## VWR Model 1235 PC

PRODUCT  
Image  
Unavailable

PRODUCT  
Image  
Unavailable

Ideal for wide variety of routine applications, designed with high-grade stainless tank in a robust outer case, features a gabled polycarbonate lid to divert condensation away from samples, over-temperature cut-out provides system protection in case bath

Rating: Not Rated Yet

**Price**

Sales price \$300.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

### Description

Ideal for wide variety of routine applications, designed with high-grade stainless tank in a robust outer case, features a gabled polycarbonate lid to divert condensation away from samples, over-temperature cut-out provides system protection in case bath accidentally operates with insufficient liquid, user adjustable secondary control system protects samples in the event of primary system failure, 11.8-Liter capacity, temperature range from ambient +5°C to 99°C, comes with base tray and gabled polycarbonate lid, 115V, 50/60Hz.