

VWR Power Supply, High Current 300 V





High current 300 V/2000 mA, 300 W maximum output.

Rating: Not Rated Yet

Price

Sales price \$625.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

Description

For high current blotting applications.

Vertical orientation to save space. High current 300 V/2000 mA, 300 W maximum output. The colour touch screen control panel is easy to use, allowing setting of voltage, amperage and run time. Four sets of output jacks allow connection of multiple electrophoresis tanks.

Specifications

Current: 2000mA

Power: 300W

Input: 100-240V

Output: 10-300V

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com