

VWR Standard Heat Block



These units are perfect for incubation and activation of cultures, enzyme reactions, and a variety of other laboratory procedures.

Rating: Not Rated Yet

Price

Sales price \$250.00

[Ask a question about this product](#)

Manufacturer [VWR](#)

Description

These units are perfect for incubation and activation of cultures, enzyme reactions, and a variety of other laboratory procedures. Interchangeable modular blocks accommodate various sample enclosures such as vials, microplates, PCR† strips, and a wide variety of tubes. Each block contains a thermometer well, and features close tube-and-block contact for maximum heat retention. Analogue temperature control knobs are easily accessible at the front of the unit and provide easy setting of minimum and maximum temperature levels.

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office

650-276-7487 fax
sales@biotechequipmentsales.com