

## Watson-Marlow 520S IP31 Peristaltic Pump



Peristaltic pump with flow-rates from 4 microliters/min up to 2.4 liters/minute; 0.1 to 220rpm.

Rating: Not Rated Yet

**Price**

Sales price \$1,850.00

[Ask a question about this product](#)

Manufacturer [Watson-Marlow, Inc](#)

### Description

Easy to install, simple to operate and inexpensive to maintain, peristaltic pumps are ideal for most fluids, including viscous, shear-sensitive, corrosive and abrasive fluids, and those containing suspended solids. Suitable for metering, dosing and dispensing applications.

520S drive offers manual control via a large-key membrane keypad: speed adjustment; run and stop; forward/reverse; "Max" key for rapid priming and purging.

Supply voltage/frequency: 100-120V/200-240V 50/60Hz

Flow rates from 4 microliters/min up to 2.4 liter/min

Precise 2,200:1 speed control range from 0.1 to 220 rpm in 0.1 rpm increments  
Large speed/flow display and easy to use tactile membrane keypad. Incorporates keypad lock  
Manual control and flow rate calibration  
MemoDose for accurate single shot dispensing

*Comes with our 90-day warranty.*

*Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.*

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office

[650-276-7487 fax](tel:650-276-7487)  
[sales@biotequipmentsales.com](mailto:sales@biotequipmentsales.com)