

## Watson-Marlow, Inc PumpPro DPM



Microprocessor control delivers rapid, precise dispensing or transfer of up to three litres per minute with a dispensing accuracy better than +/- 1%.

Rating: Not Rated Yet

**Price**

Sales price \$900.00

[Ask a question about this product](#)

Manufacturer [Watson-Marlow, Inc](#)

### Description

Microprocessor control delivers rapid, precise dispensing or transfer of up to three litres per minute with a dispensing accuracy better than +/- 1%. Membrane keypad enables operator to easily specify required volume, number of dispenses, and pump speed. Select digital display as either flowrate (mL/sec.) or speed (rpm). Dispensing may be initiated by specified time interval or with a foot switch or hand switch. Manual control available for continuous metering or transfer. Adjustable "soft-start" can be programmed to reduce splashing; end-of-dose reversing capability prevents dripping. Includes standard size, easy loading pumphead with lockable guard. Features simple calibration as well as recalibration without batch interruption.

*Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.*

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)