

## Welch 1402 vacuum pump



Welch 1402 5.6 cfm two-stage belt-drive high vacuum pump. Barely used

Rating: Not Rated Yet

**Price**

Sales price \$1,350.00

[Ask a question about this product](#)

Manufacturer [Welch](#)

Description

**!!!SOLD!!!**

This Welch 1402 two-stage belt-drive high vacuum pump is barely used. It produces 160L/min (5.6 cfm) free-air displacement and reaches  $1 \times 10^{-4}$  torr. The contoured, molded belt-guard is designed for quieter operation and eliminates visible oil weeping at the seal.