

## Yamato RE-301 Rotary Evaporator NEW IN BOX!!



\* Motorized lift \* Feedback control motor \* Voltage-sensing power supply \* Digital setting and rpm display \* Stable rotation NO BATH, NO GLASSWARE INCLUDED

Rating: Not Rated Yet

**Price**

Sales price \$1,350.00

[Ask a question about this product](#)

Manufacturer [Yamato](#)

Description	Composition		Motor	DC brushless motor (for rotation) and DC motor (for lifting)
	Performance and features	Number of revolutions		20 - 250 r/min (rpm)
		Motorized lift stroke		150mm
		Speed (rpm) setting		Control knob (digital indication)
		Safety measures (main unit)		Manual-setting lower limit, upper limit, motor overload
Lifting feature		Motorized lifting system		
Specifications	External dimension *1	With Glassware A (mm)	828(W)x 400(D)x580(H)(740(H)when raised)	
		With Glassware B/C (mm)	643(W)x 400(D)x717(H)(877(H)when raised)	
	Weight		Approx. 13kg(excluding the glassware and the water/oil bath)	
	Power supply	Main unit	100 - 240V AC, 1.5A(excluding the water/oil bath)	

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)